

—
TRAUMA



Cottey-MD

VARIABLE ANGLE
LOCKING PLATE SYSTEM

OVERFIX®
— TRAUMA —

bioadvance
— DISPOSITIVOS MÉDICOS —

anmat 
Administración Nacional de Medicamentos,
Alimentos y Tecnología Médica

CE 0123

ISO 13485

FEATURES:


- Low profile design of plate and screw minimizes the soft tissue irritation.
- Highly polished surface and rounded edges minimize the potential tendon adhesion.
- Locking and compression function can be achieved simultaneously by a locking-compression hole.



Compression




Locking



15 degree in any direction

- Innovative Variable Angle Locking Design allows both fixed angle and variable angle locking.
- Obtains better mechanical properties compared with traditional locking -compression design.
- Anatomical design makes optimal fit in bone.



- Our products are manufactured using one or more of the following materials: commercially pure titanium(ISO 5832-2, ASTM F67), Ti6Al4 V alloy (ISO 5832-3, ASTM F1472).

Cottey-MD Variable Angle Locking Plate — Mini 2.4mm

2.4 Straight Plates

Product Code	Holes	Length
52231004	4	37
52231005	5	44
52231006	6	51
52231007	7	58
52231008	8	65



2.4 L-Plates

Product Code Left	Product Code Right	Holes	Length
52231533	52231563	3	26.5
52231534	52231564	4	33.5
52231535	52231565	5	40.5
52231536	52231566	6	47.5
52231537	52231567	7	54.5



2.4 Reconstruction Plates

Product Code	Holes	Length
52231106	6	47.5
52231107	7	54.75
52231108	8	62



2.4 Condylar Plates II

Product Code Left	Product Code Right	Holes	Length
52231333	52231363	3	30
52231334	52231364	4	37.25
52231335	52231365	5	44.5
52231336	52231366	6	51.75
52231337	52231367	7	59



2.4 Condylar Plates

Product Code	Holes	Length
52231203	3	30
52231204	4	37.25
52231205	5	44.5
52231206	6	51.75
52231207	7	59



2.4 Y-Plates

Product Code	Holes	Length
52231703	3	29.7
52231704	4	36.7
52231705	5	43.7
52231706	6	50.7
52231707	7	57.7



2.4 T-Plates 3 holes head

Product Code	Holes	Length
52231603	3	27.5
52231604	4	34.5
52231605	5	41.5
52231606	6	48.5
52231607	7	55.5



2.4 T-Plates 2 holes head

Product Code	Holes	Length
52231403	3	29
52231404	4	36.25
52231405	5	43.5
52231406	6	50.75
52231407	7	58



Cottey-MD Variable Angle Locking Plate Set — Mini 2.7mm

2.7 Straight Plates

Product Code	Holes	Length
52234304	4	39
52234305	5	48
52234306	6	57
52234307	7	66
52234308	8	75
52234309	9	84
52234310	10	93



2.7 Reconstruction Plates

Product Code	Holes	Length
52234404	4	36
52234405	5	45
52234406	6	54
52234407	7	63
52234408	8	72
52234409	9	81
52234410	10	90



2.7 L-Plates

Product Code	Holes	Length
52234733	3	32
52234734	4	39.5
52234763	3	32
52234764	4	39.5



2.7 Condylar Plates

Product Code	Holes	Length
52234503	3	36
52234504	4	45
52234505	5	54
52234506	6	63
52234507	7	72
52234508	8	81



2.7 Oblique L-Plates

Product Code	Holes	Length
52234833	3	34.3
52234834	4	41.8
52234863	3	34.3
52234864	4	41.8



2.7 Condylar Plates II

Product Code Left	Product Code Right	Holes	Length
52234633	52234663	3	36
52234634	52234664	4	45
52234635	52234665	5	54
52234636	52234666	6	63
52234637	52234667	7	72
52234638	52234668	8	81



2.7 T-Plates

Product Code	Holes	Length
52234903	3	32
52234904	4	39.5



Cottey-MD Variable Angle Locking Plate Set — Wrist Solution

2.4 Distal Radius Straight Plates

Product Code	Holes	Length
52231805	5	47
52231806	6	55.5
52231807	7	64



2.4 Distal Radius Volar Plates

5 holes head

Product Code Left	Product Code Right	Holes	Length
52232033	52232063	3	43
52232035	52232065	5	61



2.4 Distal Radius Dorsal T-Plates

Product Code	Holes	Length
52231903	3	41.5
52231904	4	50.5



2.4 Distal Radius Volar Plates

8 holes head

Product Code Left	Product Code Right	Holes	Length
52232132	52232162	2	40
52232133	52232163	3	49
52232134	52232164	4	58
52232135	52232165	5	67
52232136	52232166	6	76
52232137	52232167	7	85
52232138	52232168	8	94



2.4 Distal Radius Volar extra-articular Plates 5 holes head

Product Code Left	Product Code Right	Holes	Length
52232433	52232463	3	48
52232434	52232464	4	57
52232435	52232465	5	66
52232436	52232466	6	75
52232437	52232467	7	84
52232438	52232468	8	93



2.4 Distal Radius Volar Plates 9 holes head

Product Code Left	Product Code Right	Holes	Length
52232232	52232262	2	40
52232233	52232263	3	49
52232234	52232264	4	58
52232235	52232265	5	67
52232236	52232266	6	76
52232237	52232267	7	85
52232238	52232268	8	94



2.4 Distal Radius Volar extra-articular Plates 4 holes head

Product Code Left	Product Code Right	Holes	Length
52232333	52232363	3	47
52232335	52232365	5	65



Small T-Plates (head holes 3)

Product Code	Holes	Length
52200103	3	50
52200104	4	57
52200105	5	67
52200106	6	77
52200107	7	87
52200108	8	97




2.4 Distal Radius Dorsal L-Plates 2 holes head

Product Code Left	Product Code Right	Holes	Length
52232533	52232563	3	41.5
52232534	52232564	4	50.5
52232535	52232565	5	59.5



Small T-Plates (head holes 4)

Product Code	Holes	Length
52200203	3	50
52200204	4	57
52200205	5	67
52200206	6	77
52200207	7	87
52200208	8	97



2.4 Distal Radius Dorsal L-Plates 3 holes head

Product Code Left	Product Code Right	Holes	Length
52232633	52232663	3	41.5
52232634	52232664	4	50.5
52232635	52232665	5	59.5



Oblique Radius T-Plates

Product Code Left	Product Code Right	Holes	Length
52200333	52200363	3	52
52200334	52200364	4	63
52200335	52200365	5	74
52200336	52200366	6	85
52200337	52200367	7	96
52200338	52200368	8	107



2.4 Distal Radius Dorsal Oblique L-Plates

Product Code Left	Product Code Right	Holes	Length
52232733	52232763	3	43
52232734	52232764	4	52



Cottey-MD Variable Angle Locking Plate Set — Foot Solution

2.4/2.7 I-Plates

Product Code	Holes	Length
52232802	2	17



2.4/2.7 Hind/midfoot Plates

Product Code	Holes	Length
52233105	5	31
52233108	8	45



2.4/2.7 Navicular Plates

Product Code	Holes
52232911	11



2.4/2.7 X-Plates

Product Code	Length
52233300	24
52233301	30
52233302	36



2.4/2.7 Cuboid Plates

Product Code Left	Product Code Right	Holes
52233241	52233271	11



2.4/2.7 Talus Plates

Product Code	Holes
52233006	6



2.4/2.7 Fusion Foot Plates VI

Product Code	Holes	Length
52233906	6	39
52233907	7	48



2.4/2.7 Fusion Foot Plates I

Product Code	Holes	Length
52233402	2	27
52233404	4	40



2.4/2.7 Fusion Foot Plates V

Product Code Left	Product Code Right	Holes	Length
52233832	52233862	2	44
52233833	52233863	3	53
52233834	52233864	4	62



2.4/2.7 Fusion Foot Plates II

Product Code	Length
52233500	24
52233501	30
52233502	36



2.4/2.7 Fusion Foot Plates VII-0°

Product Code Left	Product Code Right	Length
52234030	52234060	42 Small
52234031	52234061	52 Medium



2.4/2.7 Fusion Foot Plates III

Product Code	Holes	Length
52233602	2	42
52233603	3	51
52233604	4	60



2.4/2.7 Fusion Foot Plates VII-5°

Product Code Left	Product Code Right	Length
52234130	52234160	42 Small
52234131	52234161	52 Medium
52234132	52234162	57 Large



2.4/2.7 Fusion Foot Plates IV

Product Code	Holes	Length
52233702	2	46
52233703	3	55
52233704	4	64



2.4/2.7 Fusion Foot Plates VII-10°

Product Code Left	Product Code Right	Length
52234230	52234260	42 Small
52234231	52234261	52 Medium



Calcaneal Plates

Product Code Left	Product Code Right	Length
52203430	52203460	51
52203431	52203461	60
52203432	52203462	68



3.5 Calcaneal Plates II

Product Code Left	Product Code Right	Holes	Length
52205045	52205075	5	46
52205048	52205078	8	63



Cottey-MD Variable Angle Locking Plate Set — Elbow Solution

Proximal Dorsal Ulna Plates

Product Code Left	Product Code Right	Holes	Length
52200434	52200464	4	86
52200435	52200465	5	99
52200436	52200466	6	112
52200437	52200467	7	125
52200438	52200468	8	138
52200439	52200469	9	151
52200440	52200470	10	164
52200441	52200471	11	177
52200442	52200472	12	190
52200443	52200473	13	203
52200444	52200474	14	216



Distal Medial Humeral Plates

Product Code Left	Product Code Right	Holes	Length
52200633	52200663	3	58
52200634	52200664	4	71
52200635	52200665	5	84
52200636	52200666	6	97
52200637	52200667	7	110
52200638	52200668	8	123
52200639	52200669	9	136
52200640	52200670	10	149
52200641	52200671	11	162
52200642	52200672	12	175
52200643	52200673	13	188
52200644	52200674	14	201
52200645	52200675	15	214



Proximal Dorsal Ulna Plates II

Product Code	Holes	Length
52200503	3	63
52200504	4	76
52200505	5	89
52200506	6	102



Distal Posterior/Lateral Humeral Plates

Product Code Left	Product Code Right	Holes	Length
52200733	52200763	3	67
52200734	52200764	4	81
52200735	52200765	5	94
52200736	52200766	6	107
52200737	52200767	7	120
52200738	52200768	8	133
52200739	52200769	9	146
52200740	52200770	10	159
52200741	52200771	11	172
52200742	52200772	12	185
52200743	52200773	13	198
52200744	52200774	14	211



Ulna Radius Straight Plates

Product Code	Holes	Length
52203704	4	67
52203705	5	80
52203706	6	93
52203707	7	106
52203708	8	119
52203709	9	132
52203710	10	145
52203711	11	158
52203712	12	171
52203713	13	184
52203714	14	197



Cottey-MD Variable Angle Locking Plate Set — Shoulder Solution

Proximal Lateral Humeral Plates II

Product Code	Holes	Length
52200903	3	90
52200904	4	108
52200905	5	126
52200906	6	144
52200907	7	162
52200908	8	180
52200909	9	198
52200910	10	216
52200911	11	234
52200912	12	252



Superior Anterior Clavicle Plates

Product Code Left	Product Code Right	Holes	Length
52201136	52201166	6	94
52201137	52201167	7	110
52201138	52201168	8	120



Humerus Straight Plates

Product Code	Holes	Length
52203804	4	72.5
52203805	5	87.2
52203806	6	101.9
52203807	7	116.6
52203808	8	131.3
52203809	9	146
52203810	10	160.7
52203811	11	175.4
52203812	12	190.1
52203813	13	204.8
52203814	14	219.5



Superior Anterior Clavicle Plates II

Product Code Left	Product Code Right	Holes	Length
52201233	52201263	3	64
52201234	52201264	4	76
52201235	52201265	5	88
52201236	52201266	6	100
52201237	52201267	7	112
52201238	52201268	8	124
52201239	52201269	9	135



Clavicle Diaphysis Plates

Product Code Left	Product Code Right	Holes	Length
52201036	52201066	6	72
52201037	52201067	7	84
52201038	52201068	8	96
52201039	52201069	9	107.5
52201040	52201070	10	119
52201041	52201071	11	130
52201042	52201072	12	140.3



Clavicle Hook Plates

Product Code Left	Product Code Right	Holes	Length
52201334	52201364	4	15
52201335	52201365	5	15
52201336	52201366	6	15
52201337	52201367	7	15
52201344	52201374	4	18
52201345	52201375	5	18
52201346	52201376	6	18
52201347	52201377	7	18



Cottey-MD Variable Angle Locking Plate Set — Ankle Solution

1/3 Tubular Plates

Product Code	Holes	Length
52203604	4	45
52203605	5	57
52203606	6	69
52203607	7	81
52203608	8	93
52203609	9	105
52203610	10	117
52203612	12	129
52203614	14	141



3.5 Distal Medial Tibia Plates II

Product Code Left	Product Code Right	Holes	Length
52201834	52201864	4	99
52201835	52201865	5	112
52201836	52201866	6	125
52201837	52201867	7	138
52201838	52201868	8	151
52201839	52201869	9	164
52201840	52201870	10	177
52201841	52201871	11	190
52201842	52201872	12	203
52201843	52201873	13	216
52201844	52201874	14	229



Distal Lateral Fibular Plates

Product Code Left	Product Code Right	Holes	Length
52201533	52201563	3	73
52201534	52201564	4	86
52201535	52201565	5	99
52201536	52201566	6	112
52201537	52201567	7	125
52201538	52201568	8	138
52201539	52201569	9	151
52201540	52201570	10	164
52201541	52201571	11	177
52201542	52201572	12	190



3.5 Distal Lateral Tibia Plates

Product Code Left	Product Code Right	Holes	Length
52202034	52202064	4	61
52202035	52202065	5	73
52202036	52202066	6	85
52202037	52202067	7	97
52202038	52202068	8	109
52202039	52202069	9	121
52202040	52202070	10	133
52202042	52202072	12	157
52202044	52202074	14	181
52202046	52202076	16	205
52202048	52202078	18	229
52202050	52202080	20	253



3.5 Distal Medial Tibia Plates

Product Code Left	Product Code Right	Holes	Length
52201734	52201764	4	109
52201735	52201765	5	121
52201736	52201766	6	133
52201737	52201767	7	145
52201738	52201768	8	157
52201739	52201769	9	169
52201740	52201770	10	181
52201741	52201771	11	193
52201742	52201772	12	205
52201743	52201773	13	217
52201744	52201774	14	229



4.5 Proximal Medial Tibia Plates

Product Code Left	Product Code Right	Holes	Length
52202534	52202564	4	81
52202535	52202565	5	97
52202536	52202566	6	113
52202537	52202567	7	129
52202538	52202568	8	145
52202539	52202569	9	161
52202540	52202570	10	177



3.5 Metaphyseal Plates

Product Code	Holes	Length
52204206	86	6
52204207	99	7
52204208	112	8
52204209	125	9
52204210	138	10
52204211	151	11
52204212	164	12
52204214	190	14
52204216	216	16
52204218	242	18



3.5 Reconstruction Plates

Product Code	Holes	Length
52204404	4	56
52204405	5	70
52204406	6	84
52204407	7	98
52204408	8	112
52204409	9	126
52204410	10	140
52204412	12	168
52204414	14	196
52204416	16	224



Cotley-MD Variable Angle Locking Plate Set — Knee Solution

4.5/3.5 Metaphyseal Plates

Product Code	Holes	Length
52204308	8	119
52204309	9	136
52204310	10	153
52204311	11	170
52204312	12	187
52204313	13	204
52204314	14	221
52204316	16	255
52204318	18	289
52204320	20	323



4.5 Proximal Lateral Tibia Plates

Product Code Left	Product Code Right	Holes	Length
52202334	52202364	4	82
52202335	52202365	5	100
52202336	52202366	6	118
52202337	52202367	7	136
52202338	52202368	8	154
52202339	52202369	9	172
52202340	52202370	10	190
52202341	52202371	11	208
52202342	52202372	12	226
52202344	52202374	14	262
52202346	52202376	16	298



Femur Curved Plates

Product Code	Holes	Length
52204106	6	143
52204107	7	161
52204108	8	179
52204109	9	197
52204110	10	215
52204111	11	233
52204112	12	251
52204113	13	269
52204114	14	287
52204115	15	305
52204116	16	323
52204117	17	341
52204118	18	359



Distal Lateral Femur Plates

Product Code Left	Product Code Right	Holes	Length
52202736	52202766	6	158
52202737	52202767	7	176
52202738	52202768	8	194
52202739	52202769	9	212
52202740	52202770	10	230
52202741	52202771	11	248
52202742	52202772	12	266
52202743	52202773	13	284
52202744	52202774	14	302



3.5/4.5 Proximal Lateral Tibia Plates

Product Code Left	Product Code Right	Holes	Length
52202134	52202164	4	93
52202135	52202165	5	111
52202136	52202166	6	129
52202137	52202167	7	147
52202138	52202168	8	165
52202139	52202169	9	183
52202140	52202170	10	201
52202141	52202171	11	219
52202142	52202172	12	237
52202144	52202174	14	273
52202146	52202176	16	309



3.5 Proximal Lateral Tibia Plates

Product Code Left	Product Code Right	Holes	Length
52202234	52202264	4	81
52202235	52202265	5	94
52202236	52202266	6	107
52202237	52202267	7	120
52202238	52202268	8	133
52202239	52202269	9	146
52202240	52202270	10	159
52202241	52202271	11	172
52202242	52202272	12	185
52202244	52202274	14	211
52202246	52202276	16	237



Tibia Straight Plates

Product Code	Holes	Length
52203904	4	85.5
52203905	5	103
52203906	6	120.5
52203907	7	138
52203908	8	155.5
52203909	9	173
52203910	10	190.5
52203911	11	208
52203912	12	225.5
52203913	13	243
52203914	14	278
52203916	16	313
52203918	18	348



Cotley-MD Variable Angle Locking Plate Set — Hip Solution

Proximal Lateral Femur Plates

Product Code Left	Product Code Right	Holes	Length
52202834	52202864	4	118
52202835	52202865	5	136
52202836	52202866	6	154
52202837	52202867	7	172
52202838	52202868	8	190
52202839	52202869	9	208
52202840	52202870	10	226
52202841	52202871	11	244
52202842	52202872	12	262
52202843	52202873	13	280



Femur Straight Plates

Product Code	Holes	Length
52204006	6	133
52204007	7	151
52204008	8	169
52204009	9	187
52204010	10	205
52204011	11	223
52204012	12	241
52204013	13	259
52204014	14	277
52204016	16	313
52204018	18	349



Cottey-MD Variable Angle Locking Plate Set — Screw

2.4mm Variable Angle Locking Screw

Product Code	Length
53002406	6
53002407	7
53002408	8
53002409	9
53002410	10
53002411	11
53002412	12
53002413	13
53002414	14
53002416	16
53002418	18
53002420	20
53002422	22
53002424	24
53002426	26
53002428	28
53002430	30



2.7mm Variable Angle Locking Screw

Product Code	Length
53002708	8
53002710	10
53002712	12
53002714	14
53002716	16
53002718	18
53002720	20
53002722	22
53002724	24
53002726	26
53002728	28
53002730	30
53002732	32
53002734	34
53002736	36
53002738	38
53002740	40
53002742	42
53002744	44
53002746	46
53002748	48
53002750	50
53002755	55
53002760	60



2.4mm Cortex Screw

Product Code	Length
52402406	6
52402408	8
52402410	10
52402412	12
52402414	14
52402416	16
52402418	18
52402420	20
52402422	22
52402424	24
52402426	26
52402428	28
52402430	30
52402432	32
52402434	34
52402436	36
52402438	38
52402440	40
52402442	42
52402444	44
52402446	46
52402448	48
52402450	50



2.7 Locking Screw

Product Code	Length
52602708	8
52602710	10
52602712	12
52602714	14
52602716	16
52602718	18
52602720	20
52602722	22
52602724	24
52602726	26
52602728	28
52602730	30
52602732	32
52602734	34
52602736	36
52602738	38
52602740	40
52602742	42
52602744	44
52602746	46
52602748	48
52602750	50
52602755	55
52602760	60



2.7 Cortex Screw

Product Code	Length
52402706	6
52402708	8
52402710	10
52402712	12
52402714	14
52402716	16
52402718	18
52402720	20
52402722	22
52402724	24
52402726	26
52402728	28
52402730	30
52402732	32
52402734	34
52402736	36
52402738	38
52402740	40
52402742	42
52402744	44
52402746	46
52402748	48
52402750	50
52402755	55
52402760	60



3.5 Cortex Screw

Product Code	Length
52403510	10
52403512	12
52403514	14
52403516	16
52403518	18
52403520	20
52403522	22
52403524	24
52403526	26
52403528	28
52403530	30
52403532	32
52403534	34
52403536	36
52403538	38
52403540	40
52403542	42
52403544	44
52403546	46
52403548	48
52403550	50
52403555	55
52403560	60
52403565	65
52403570	70
52403575	75
52403580	80



3.5mm Variable Angle Locking Screw

Product Code	Length
53003510	10
53003512	12
53003514	14
53003516	16
53003518	18
53003520	20
53003522	22
53003524	24
53003526	26
53003528	28
53003530	30
53003532	32
53003534	34
53003536	36
53003538	38
53003540	40
53003542	42
53003544	44
53003546	46
53003548	48
53003550	50
53003555	55
53003560	60
53003565	65
53003570	70
53003575	75
53003580	80
53003585	85
53003590	90
53003595	95



5.0mm Variable Angle Locking Screw

Product Code	Length
53005014	14
53005016	16
53005018	18
53005020	20
53005022	22
53005024	24
53005026	26
53005028	28
53005030	30
53005032	32
53005034	34
53005036	36
53005038	38
53005040	40
53005042	42
53005044	44
53005046	46
53005048	48
53005050	50
53005055	55
53005060	60
53005065	65
53005070	70
53005075	75
53005080	80
53005085	85
53005090	90



5.5 Cancellous Screw Fully Threaded

Product Code	Length
52505530	30
52505535	35
52505540	40
52505545	45
52505550	50
52505555	55
52505560	60
52505565	65
52505570	70
52505575	75
52505580	80
52505585	85
52505590	90
52505595	95
52505500	100
52505505	105
52505510	110



5.5 Cancellous Screw

Product Code	Length	Thread length	Product Code	Length	Thread length	Product Code	Length	Thread length
52525530	30	16	52535540	40	24	52545550	50	32
52525535	35	16	52535545	45	24	52545555	55	32
52525540	40	16	52535550	50	24	52545560	60	32
52525545	45	16	52535555	55	24	52545565	65	32
52525550	50	16	52535560	60	24	52545570	70	32
52525555	55	16	52535565	65	24	52545575	75	32
52525560	60	16	52535570	70	24	52545580	80	32
52525565	65	16	52535575	75	24	52545585	85	32
52525570	70	16	52535580	80	24	52545590	90	32
52525575	75	16	52535585	85	24	52545595	95	32
52525580	80	16	52535590	90	24	52545500	100	32
52525585	85	16	52535595	95	24	52545505	105	32
52525590	90	16	52535500	100	24	52545510	110	32
52525595	95	16	52535505	105	24			
52525500	100	16	52535510	110	24			
52525505	105	16						
52525510	110	16						



4.0mm Cancellous Screws

Product Code	Length	Thread length
52514010	10	5
52514012	12	5
52514014	14	5
52514016	16	6
52514018	18	7
52514020	20	8
52514022	22	9
52514024	24	10
52514026	26	12
52514028	28	14
52514030	30	14
52514032	32	14
52514034	34	14
52514036	36	14
52514038	38	14
52514040	40	14
52514045	45	14
52514050	50	14
52514055	55	14
52514060	60	14
52514065	65	14
52514070	70	14
52514075	75	14
52514080	80	14
52514085	85	14
52514090	90	14
52514095	95	14
52514000	100	14



4.0mm Cancellous Screws, Fully Threaded

Product Code	Length
52504010	10
52504012	12
52504014	14
52504016	16
52504018	18
52504020	20
52504022	22
52504024	24
52504026	26
52504028	28
52504030	30
52504032	32
52504034	34
52504036	36
52504038	38
52504040	40
52504045	45
52504050	50
52504055	55
52504060	60
52504065	65
52504070	70
52504075	75
52504080	80
52504085	85
52504090	90
52504095	95
52504000	100



4.5mm Cortex Screw

Product Code	Length
52404514	14
52404516	16
52404518	18
52404520	20
52404522	22
52404524	24
52404526	26
52404528	28
52404530	30
52404532	32
52404534	34
52404536	36
52404538	38
52404540	40
52404542	42
52404544	44
52404546	46
52404548	48
52404550	50
52404552	52
52404554	54
52404556	56
52404558	58
52404560	60
52404562	62
52404564	64
52404566	66
52404568	68
52404570	70
52404575	75
52404580	80
52404585	85
52404590	90



6.5mm Locking Screws

Product Code	Length
52606530	30
52606535	35
52606540	40
52606545	45
52606550	50
52606555	55
52606560	60
52606565	65
52606570	70
52606575	75
52606580	80
52606585	85
52606590	90
52606595	95
52606500	100
52606505	105
52606510	110
52606515	115
52606520	120
52606525	125



5.5mm Variable Angle Locking Screw Fully Threaded

Product Code	Length
53005530	30
53005535	35
53005540	40
53005545	45
53005550	50
53005555	55
53005560	60
53005565	65
53005570	70
53005575	75
53005580	80
53005585	85
53005590	90
53005595	95
53005500	100
53005505	105
53005510	110
53005515	115
53005520	120
53005525	125



5.5mm Variable Angle Locking Screw

Product Code	Length	Thread length
53015530	30	16
53015535	35	16
53015540	40	16
53015545	45	16
53015550	50	16
53015555	55	24
53015560	60	24
53015565	65	24
53015570	70	24
53015575	75	24
53015580	80	32
53015585	85	32
53015590	90	32
53015595	95	32
53015500	100	32
53015505	105	32
53015510	110	32
53015515	115	32
53015520	120	32
53015525	125	32




Desde nuestra sede operativa en Estados Unidos, y filiales en México, Argentina, Bolivia y Chile, ofrecemos las gamas más completas para trauma, columna, cadera y rodilla.



Nuestros dispositivos cuentan con registros sanitarios y cumplen las certificaciones de calidad vigentes en cada país. Implementamos procedimientos de trazabilidad y tecnovigilancia en las diferentes etapas de almacenamiento y comercialización. Esta rigurosidad en el control nos permite sostener un alto valor de marca, garantizando la mayor seguridad a los pacientes.

Con más de 8 años dedicados a la importación y comercialización de dispositivos médicos, nos destacamos por nuestro servicio, productos de eficiencia clínica y costos competitivos. Para conocer nuestras líneas de productos, por favor póngase en contacto para vincularlo con el distribuidor más cercano a su ubicación.



USA
1001 N Federal Hwy
S. 355, Hallandale Beach
FL 33009
+1786 375 3968
bioadvanceusa.com

CDMX
Amores 1322
Colonia del Valle Centro
CP 03100, CDMX
+52 55 5925 5323
bioadvance.com.mx

GUADALAJARA
Av. Patria 179. Of. 302 y 303, Piso 3
CP. 45030. Col. Prados Guadalupe
Zapopan, Jalisco
+52 33 1656 4268
bioadvance.com.mx

ARGENTINA
Cnel. Cetz 336, Piso 2
Martínez, B1642
Buenos Aires
+54 11 4765 5995
bioadvance.com.ar

BOLIVIA
Calle Pedro Rodríguez 3830
Santa Cruz de la Sierra
Bolivia
+59 1 7048 5050
bolivia.bioadvanceusa.com

CHILE
Av. Alonso de Cordova 5870
Of. 1711, 7560885
Santiago de Chile
+56 9 8493 9874
chile.bioadvanceusa.com



OVERFIX[®]

— TRAUMA —

bioadvance

— DISPOSITIVOS MÉDICOS —

MÉXICO

bioadvance.com.mx

ARGENTINA

bioadvance.com.ar

BOLIVIA

bolivia.bioadvanceusa.com

CHILE

chile.bioadvanceusa.com